

PROAIR® DIGIHALER™ INHALER

DRY-POWDER INHALER

How to Use Your ProAir® Digihaler™ Inhaler

Your clinician has prescribed medication that comes in a ProAir Digihaler inhaler. It is important that you know how to use the ProAir Digihaler device so you can take your medication correctly and breathe better.

Important: The dose indicator window on the back of the inhaler shows the number of doses left in the inhaler. The device counts down by 1 each time you open the red cap.

Note: You do not need to shake or prime this inhaler.


To use your ProAir® Digihaler™ Inhaler, complete these steps:




1. With the red cap closed, hold the inhaler upright.



2. Open the red cap. You should hear a click.



3. Remove any dust or lint from the mouthpiece using a clean cloth or tissue.
Important: Do not wash or put the inhaler in water.



4. Hold the inhaler away from your face, and then breathe out completely.
Important: Do not blow into the inhaler.




5. Hold the inhaler upright, place the mouthpiece between your lips, and then close your lips around the mouthpiece.



6. Take one long, steady, deep breath through your mouth and the inhaler. Breathe in fast and deep. **Important:** Do not block the vents above the mouthpiece with your lips or fingers.



7. Remove the inhaler from between your lips.
8. Hold your breath for 10 seconds (or as long as you can).



9. Turn your head, and breathe out slowly and gently away from the inhaler.



10. Close the red cap.

11. If your clinician has told you to take more than one dose of medication, repeat steps 2-9.

When fewer than 20 doses remain in the ProAir Digihaler device, the counter in the dose indicator window turns red to remind you to get a refill. When you have inhaled the last dose in the ProAir Digihaler device, the dose indicator window will show 0. Your inhaler is empty. Dispose of the inhaler properly. Keep the ProAir Digihaler device and all medications out of the reach of children and pets.

