

How to Use Your Small Volume Compressor Nebulizer

To make your breathing better, you **MUST** take your medicine as explained below. Following these instructions puts more of the medicine into your lungs and help keep your nebulizer clean. This will open up your air passages and help you breathe easier and feel better.

To use your small volume compressor nebulizer, complete these steps:



1. Place compressor on a level steady surface, plug power cord into a wall outlet.



2. Connect one end of the nebulizer tubing to the air outlet connector.



3. Wash hands and add prescribed medication to the nebulizer using dropper or unit dose vial.



4. Assemble mouthpiece, T-piece, and reservoir tubing to the nebulizer.






5. Attach tubing to nebulizer air inlet.



6. Press power switch to "on."



7. Sit upright in a comfortable position. Place mouthpiece between teeth, over the tongue, with lips comfortably sealed. Take deep relaxed breaths through your mouth until nebulizer begins to sputter.



8. Turn off power to compressor.