

PROAIR RESPICLICK® INHALER

DRY-POWDER INHALER

How to Use Your ProAir RespiClick® Inhaler

Your clinician has prescribed medication that comes in a RespiClick inhaler. It is important that you know how to use the RespiClick device so that you can take your medication correctly and breathe better.

Important: Before using your RespiClick inhaler for the first time, check the dose indicator window on the back of the device to determine how many doses are in your inhaler. The inhaler counts down by 1 each time you use the inhaler.

To use your ProAir RespiClick® Inhaler, complete these steps:



1. Hold the inhaler upright with the clear mouthpiece facing you.



2. Pull the clear mouthpiece cover back until you hear a click.



3. Hold the inhaler away from your face, and then breathe out. **Important:** Do not blow into the inhaler.




4. Place the mouthpiece between your lips, then close your lips around the mouthpiece. **Important:** Do not block the air vent above the mouthpiece with your lips or fingers.
5. Breathe in quickly and deeply through your mouth and the inhaler.



6. Remove the inhaler from between your lips.
7. Hold your breath for 10 seconds. If you cannot hold your breath for 10 seconds, hold it as long as you can.



8. Turn your head away from the inhaler and breathe out slowly.



9. Check the dose indicator window on the back of your inhaler to make sure a dose was delivered.



10. Firmly close the clear cap over the mouthpiece.
11. If your clinician has told you to take more than one dose, repeat steps 2-8.

- ! Store your RespiClick inhaler in a cool, dry place. To clean your inhaler, gently wipe it with a clean cloth or tissue.
- Keep your RespiClick inhaler and all medications out of the reach of children and pets.

